


















ShapeWorks Workflow

Section 1: Installation

- To get ShapeWorks, first [install miniconda](#) if you don't already have it.
- Visit ShapeWorks's [GitHub](#) and download the latest stable release.
- Install ShapeWorks using the installer you downloaded
- Open a miniconda terminal and navigate to the ShapeWorks folder
 - Default location is C:\Program Files\ShapeWorks
- Run install_shapeworks in the miniconda terminal
 - NOTE: If this step takes a while to complete, try deleting the shapeworks folder and updating your miniconda
- Once that is complete, navigate to C:\Program Files\ShapeWorks\bin and run ShapeWorksStudio.exe

Section 2: Importing a dataset

- For our purposes, we are using ShapeWorks to analyze .stl files.
- Download our shoulder dataset from [here](#)
- Once it is downloaded, clone the nick-dev branch from the Shoulder-Keypoint repo on GitHub and create a folder named "Projects" in the root of that folder

 .git	2/25/2022 2:43 PM	File folder	
 config	2/15/2022 8:32 PM	File folder	
 Projects	2/25/2022 3:24 PM	File folder	
 .gitignore	2/16/2022 12:41 AM	GITIGNORE File	1 KB
 ____utility.py	2/15/2022 8:32 PM	Python File	7 KB
 ____etl.py	2/15/2022 8:32 PM	Python File	3 KB
 ____train.py	2/15/2022 8:32 PM	Python File	9 KB
 ____test.py	2/15/2022 8:32 PM	Python File	9 KB
 CreateProject.py	2/25/2022 3:24 PM	Python File	5 KB
 FeaturePointDataset.py	2/25/2022 2:11 PM	Python File	26 KB
 FeaturePointDataset_tester.py	2/25/2022 2:11 PM	Python File	22 KB
 jensen_notes_11.4.21.txt	2/15/2022 8:32 PM	TXT File	1 KB
 PointHRNet.ipynb	2/15/2022 8:32 PM	IPYNB File	64 KB
 pose_hrnet.py	2/15/2022 8:32 PM	Python File	20 KB
 README.md	2/15/2022 8:32 PM	MD File	1 KB
 requirements.txt	2/15/2022 8:32 PM	TXT File	1 KB
 tester.ipynb	2/25/2022 2:11 PM	IPYNB File	0 KB

- Open the repo in VSCode and run the file "CreateProject.py"
 - The program will prompt you first for a name for the project file.
 - Next, get the path to the shoulder data you downloaded earlier and paste it into the program. It should be that path to a folder that contains the two datasets, "Akira_Organized" and "Keisuke_Organized"
 - Next, choose which bones to import. Type 1 for scapula, 2 for humerus, or 3 for clavicle.
 - Lastly, choose the amount of bones to import. I would recommend 9 at most if you're just trying to get optimization to work correctly and mess around with parameters.
 - If you're trying to generate data for the whole dataset, however, type "all"
 - The project file will be saved in the Projects folder you created earlier
- Open ShapeWorks and click "Open Existing Project", then navigate to the Projects folder and select the file you created above to import your desired .stl files.

Section 3: Using ShapeWorks to create point clouds

- ShapeWorks uses three steps in order to generate a point cloud:
 - *Groom*, which smooths out the models and makes them easier for the program to process
 - *Optimize*, which is what we use to actually generate the points on each mesh
 - *Analyze*, which we can use to visualize the different point clouds for each individual mesh and see how well our optimization parameters did

Section 3a: Grooming


Mesh Grooming

- Mesh Grooming has three main options:
 - Fill Holes
 - Smooth
 - Remesh
- Fill holes should always be selected for our purposes, as many of the scapula meshes have...you guessed it: holes!
 - NOTE: If the grooming step is taking a long time for Scapulae, de select this option. Not recommended, but it was not completing on some computers so do this if needed.
- Smoothing is pretty self explanatory: it smooths out the mesh. This *can* be useful, but if you enable it be sure to keep iteration amounts between 1 and 3, as going higher loses the surface detail that we need to generate accurate clouds.
- Remesh modifies the mesh so that the triangles that make it up have more similar areas. This is very useful, but keep it at default.

Alignment

- This section allows you to modify the alignment and the reflections of meshes
- For our purposes, the *CreateProject.py* file mentioned above assigns the "right" group to any bone that contains "Raxes" in its name. Otherwise, it assigns the "left" group
 - This is done because the Keisuke dataset is the only one with meshes from both sides of the body - Akira has only the left. The Keisuke dataset names corresponding meshes "Raxes" for the right side, and "Laxes" for the left, resulting in the above logic for the script.
- Select the "Reflect" option to reflect the meshes so that they all have the same orientation.

- Change the "shape_file" option to "group_side" to enable this feature




Section 3b: Optimizing

- There are a lot of options to choose from here, and they all have different effects. Here is what I have found through trial and error, and [the page on what the parameters mean](#) and some [quick tips on optimization](#) from ShapeWorks:
 - Number of particles
 - The amount of particles the program generates
 - Keep this number low while fine-tuning your parameters ($x < 64$) but increase and decrease by factors of $x \times 2$ if you need to.
 - Increasing it provides a mesh with higher detail at the end, but takes much longer
 - Initial Relative Weighting
 - How the initial correspondence is weighted
 - Increase this number (by a factor of 0.01) in order to make particles correspond better
 - Relative Weighting
 - Relative weighting of correspondence during optimization
 - Decrease this number (by a factor of 1) to help particles spread more evenly across the mesh
 - Starting Regularization & Ending Regularization
 - Don't worry about modifying these(?)
 - Iterations Per Split
 - How many times the optimization occurs per split of points
 - Keep this around 1000, increase it to help get particles to be in the "right places" but there are diminishing returns on increasing it.
 - This affects processing time quite a lot, especially with higher amounts of particles
 - Optimization Iterations
 - How many times the program will optimize point location per iteration. Keep this around 1000.
 - Geodesic Distance

- Don't enable this. Increases computing time by a factor of 10x and is not useful for our application.
- Normals
 - Tells the algorithm whether to consider surface normals and the xyz of a particle in correspondence
 - Should be enabled to avoid particles flipping sides on thin structures
 - Seems like it may be needed in order to get particles to place correctly on most smaller structures
- Normals Strength
 - How much the algorithm ascribes strength of surface normals relative to position.
- Procrustes
 - These options provide minimal changes and have little documentation. Would recommend to just avoid for now.
- Multiscale Mode
 - No documentation on this option and it hasn't provided meaningful improvements. Would recommend avoiding this.
- Multiscale Start
 - After how many particles should multiscale mode begin.
- Narrow Band
 - Doesn't affect optimization. Don't worry about this.

Section 3c: Exporting Data

- To export the key point data, just save the project! This will save all of the current particles in a folder with the same name as the project, just with "_particles" appended.



Name	Date modified	Type	Size
all-clav_particles	3/22/2022 3:21 PM	File folder	
all-hums_particles	3/22/2022 3:21 PM	File folder	
all-scaps_particles	3/17/2022 10:35 AM	File folder	
groomed	3/22/2022 3:16 PM	File folder	
all-clav.xlsx	3/22/2022 3:54 PM	Microsoft Excel W...	6 KB
all-hums.xlsx	3/22/2022 3:54 PM	Microsoft Excel W...	6 KB
all-scaps.xlsx	3/22/2022 3:54 PM	Microsoft Excel W...	6 KB